

Product Description SALSA[®] MLPA[®] Probemix P214-C2 COL2A1

To be used with the MLPA General Protocol.

Version C2

As compared to version C1, one reference probe has been replaced, one reference probe has been removed, and the sequence of one probe has been adjusted. For complete product history see page 8.

Catalogue numbers:

- P214-025R: SALSA MLPA Probemix P214 COL2A1, 25 reactions.
- P214-050R: SALSA MLPA Probemix P214 COL2A1, 50 reactions.
- P214-100R: SALSA MLPA Probemix P214 COL2A1, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

The SALSA MLPA Probemix P214 COL2A1 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *COL2A1* gene, which is associated with type II collagen disorders.

The *COL2A1* gene (54 exons) spans ~31.5 kb of genomic DNA and is located on 12q13.11, ~47 Mb from the p-telomere. It encodes the pro-alpha-1 (II) chain of type II collagen, the major collagen found in cartilage and the vitreous humour of the eye. Defects in the *COL2A1* gene lead to a number of different disorders, collectively called type II collagen disorders, which include amongst others achondrogenesis type 2, Stickler syndrome type 1, Kniest dysplasia, and spondyloepiphyseal dysplasia congenita. These disorders are characterised by abnormalities in the ocular, skeletal, orofacial, and audiological systems.

Mutations in *COL2A1* have also been identified in individuals with chondrosarcoma, the second most common primary bone malignancy after osteosarcoma. Tarpey et al. (2013) reported insertions, deletions and rearrangements in 37% of the chondrosarcoma cases.

More information is available at https://www.ncbi.nlm.nih.gov/books/NBK540447/ (type II collagen disorders) and https://www.ncbi.nlm.nih.gov/books/NBK540447/ (type II collagen disorders) and https://www.ncbi.nlm.nih.gov/books/NBK540447/ (type II collagen disorders) and https://www.ncbi.nlm.nih.gov/books/NBK1302/ (Stickler syndrome).

This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene For NM_ mRNA reference sequences: http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide Locus Reference Genomic (LRG) database: http://www.lrg-sequence.org/

Exon numbering

The COL2A1 exon numbering used in this P214-C2 COL2A1 product description is the exon numbering from the NG_008072.1 sequence. The exon numbering of the NM_ sequence that was used for determining a

probe's ligation site does not always correspond to the exon numbering obtained from the NG sequences. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P214-C2 COL2A1 contains 57 MLPA probes with amplification products between 124 and 500 nucleotides (nt). This includes 46 probes covering 43 different exons out of 54 exons of the *COL2A1* gene. In addition, 11 reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name			
64-70-76-82	64-70-76-82 Q-fragments (only visible with <100 ng sample DNA)			
88-96	D-fragments (low signal indicates incomplete denaturation)			
92	Benchmark fragment			
100	X-fragment (X chromosome specific)			
105	Y-fragment (Y chromosome specific)			

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of type II collagen disorders. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (https://catalog.coriell.org) and Leibniz Institute DSMZ (https://www.dsmz.de/) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.



Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	0.80 < FR < 1.20
Homozygous deletion	FR = 0
Heterozygous deletion	0.40 < FR < 0.65
Heterozygous duplication	1.30 < FR < 1.65
Heterozygous triplication/homozygous duplication	1.75 < FR < 2.15
Ambiguous copy number	All other values

Note: The term "dosage quotient", used in older product description versions, has been replaced by "final ratio" to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- <u>Arranging probes</u> according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- <u>False positive results</u>: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- <u>Normal copy number variation</u> in healthy individuals is described in the database of genomic variants: <u>http://dgv.tcag.ca/dgv/app/home</u>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- <u>Copy number changes detected by reference probes</u> or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more
 exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale
 peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net
 software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun



the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure

- In most populations, the major cause of genetic defects in the *COL2A1* gene are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P214 COL2A1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false
 positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe
 signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe
 oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

COL2A1 mutation database

https://databases.lovd.nl/shared/genes/COL2A1. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on http://varnomen.hgvs.org/.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *COL2A1* exons 2 and 4 but not exon 3) to MRC Holland: info@mrcholland.com.



Table 1. SALSA MLPA Probemix P214-C2 COL2A1

ength (nt)	SALSA MLPA probe		position (hg18)ª
	-	Reference	COL2A1
64-105	Control fragments – see table in probemix	content section for more info	rmation
124	Reference probe 19616-L26241	4p	
130	COL2A1 probe 20226-L27557		Exon 40
139	COL2A1 probe 07387-L27860		Exon 1
145	COL2A1 probe 07397-L27861		Exon 24
150	COL2A1 probe 20227-L28154		Exon 33
154	COL2A1 probe 07393-L29025		Exon 16
160	Reference probe 02417-L27862	бр	
165	COL2A1 probe 20228-L28265		Exon 26
170	COL2A1 probe 07394-L28264		Exon 17
175	COL2A1 probe 15449-L28263		Exon 10
179 *	Reference probe 01779-L17334	13q	
184	COL2A1 probe 20229-L27560		Exon 18
190	COL2A1 probe 20230-L27561		Exon 6
196	Reference probe 18769-L27864	Зр	
201	COL2A1 probe 15450-L27865		Exon 23
207	COL2A1 probe 20232-L27866		Exon 12
213	COL2A1 probe 07404-L27958		Exon 46
218	COL2A1 probe 15291-L27868		Exon 1
225	COL2A1 probe 07405-L28736		Exon 49
229	COL2A1 probe 15451-L27870		Exon 14
234	Reference probe 16398-L27959	17q	
238	COL2A1 probe 07406-L27872	1	Exon 51
244	COL2A1 probe 20233-L27564		Exon 4
252	COL2A1 probe 15452-L25045		Exon 25
257	COL2A1 probe 15453-L25046		Exon 29
265 Ж	COL2A1 probe 20234-SP0942-L27565		Exon 22
273 Δ	COL2A1 probe 20235-L27566		Exon 37
278	Reference probe 18594-L28268	2q	Exer ex
283	COL2A1 probe 07390-L27874	24	Exon 8
289	COL2A1 probe 07390 127074		Exon 35
298	COL2A1 probe 13289-L27873		Exon 27
298 305 Ж	COL2A1 probe 07398-L27877 COL2A1 probe 20463-SP0949-L28266		Exon 50
	COL2A1 probe 20403-3P0949-L28200 COL2A1 probe 15455-L27878		
312	Reference probe 08987-L28685	0.7	Exon 53
319		9q	Even EA
325	COL2A1 probe 07407-L28267		Exon 54
331	COL2A1 probe 20238-L27569		Exon 7
<u>337 Ж</u>	COL2A1 probe 20239-SP0944-L27954		Exon 34
346 Ж	COL2A1 probe 20241-SP0946-L27572		Exon 48
353	COL2A1 probe 07402-L27881		Exon 39
361 Ж	COL2A1 probe 20242-SP0947-L27882		Exon 41
368	COL2A1 probe 07863-L27883		Exon 13
374	COL2A1 probe 15456-L27884	-	Exon 52
382	Reference probe 17429-L27885	8p	
395 Ж	COL2A1 probe 20243-SP0948-L28915		Exon 44
403	COL2A1 probe 20244-L28689		Exon 9
409	COL2A1 probe 20245-L27576		Exon 42
418	COL2A1 probe 15459-L28687		Exon 31
427	COL2A1 probe 15460-L28026		Exon 54
434	Reference probe 08839-L27962	2р	
441	COL2A1 probe 20246-L27963		Exon 20
445	COL2A1 probe 20247-L27964		Exon 43



Longth (nt)		Chromosomal position (hg18) ^a		
Length (nt)	SALSA MLPA probe	Reference	COL2A1	
454	COL2A1 probe 20248-L27579		Exon 11	
463	COL2A1 probe 20249-L27580		Exon 2	
476	Reference probe 20551-L28686	11q		
485 ¥	COL2A1 probe 20252-L32985		Exon 3	
494	COL2A1 probe 20253-L27584		Exon 50	
500	Reference probe 14894-L27890	15q		

^a See section Exon numbering on page 1 for more information.

* New in version C2.

¥ Changed in version C2. Minor alteration, no change in sequence detected.

 Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

X This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 2. COL2A1	probes arrang	ged according	to chromosomal location
-----------------	---------------	---------------	-------------------------

Length (nt)	SALSA MLPA probe	COL2A1 exonª	Ligation site NM_001844.5	<u>Partial</u> sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	156-158 (Exon 1)		
218	15291-L27868	Exon 1	387 nt before exon 1	CGTCCTTGGTCT-AGGGCTCTCGGC	0.7 kb
139	07387-L27860	Exon 1	40 nt after exon 1	GCCTGCTTTCCA-TGCGTCCCTCAG	4.3 kb
463	20249-L27580	Exon 2	444-445	TCGCCACTGCCA-GTGGTTGTAATT	1.5 kb
485	20252-L32985	Exon 3	455-456	GCAGGGCAACCA-GGACCAAAGGTA	0.2 kb
244	20233-L27564	Exon 4	494-495	GGAGACATCAAG-GATGTAAGTGCA	0.3 kb
190	20230-L27561	Exon 6	17 nt after exon 6, reverse	AGGGTCAAGCAG-CATTGCTTTTA	0.3 kb
331	20238-L27569	Exon 7	19 nt after exon 7, reverse	GCCTGAAGGAAT-GGGAAGTAAGGA	1.0 kb
283	07390-L27874	Exon 8	722-723	GGATTTGATGAA-AAGGCTGGTGGC	0.7 kb
403	20244-L28689	Exon 9	802-801, reverse	TTACAGGAGCAC-CTGCAGGGCCTG	0.1 kb
175	15449-L28263	Exon 10	13 nt before exon 10	CTGGTATCCTCA-TTTTACTTTTA	0.5 kb
454	20248-L27579	Exon 11	913-914	AAAGCCTGGTGA-TGATGTGAGTAT	0.8 kb
207	20232-L27866	Exon 12	30 nt before exon 12	GCGTCTCTGAGG-AAGCTGGGATAT	0.5 kb
368	07863-L27883	Exon 13	992-993	GGTTTCCCAGGA-ACCCCAGGCCTT	0.2 kb
229	15451-L27870	Exon 14	2 nt before exon 14	CTCTTGTTCCCT-AGGGTTATCCAG	1.0 kb
154	07393-L29025	Exon 16	1147-1148	CCTGCCTGGTGA-AAGAGGACGGAC	3.1 kb
170	07394-L28264	Exon 17	9 nt before exon 17	TCACTTCCTTCT-TGCTCACAGGGT	0.6 kb
184	20229-L27560	Exon 18	1270-1271	TGGTGCTCCTGG-AGCCAAGGTACG	1.9 kb
441	20246-L27963	Exon 20	1397-1398	GGAACAGATGGA-ATTCCTGGAGCC	0.5 kb
265 Ж	20234-SP0942- L27565	Exon 22	1550-1551 and 6 nt after exon 22	GGCTTCAAAGGT-30 nt spanning oligo-ATCTGCCCCCAA	0.5 kb
201	15450-L27865	Exon 23	9 nt after exon 23	AGAGTTAAGTGA-ATGTGGAGGCTC	0.4 kb
145	07397-L27861	Exon 24	1735-1736	GGCAGGTCCCAA-GGTGAGTGGGAG	0.1 kb
252	15452-L25045	Exon 25	3 nt before exon 25	TGTGTACCCTTG-TAGGGAGCCCCT	0.3 kb
165	20228-L28265	Exon 26	10 nt after exon 26, reverse	AGCCCTCAGAGG-ATAGACTTACAG	0.5 kb
298	07398-L27877	Exon 27	1919-1920	GGTCGTCCTGGA-CCTCCAGGTCCT	0.8 kb
257	15453-L25046	Exon 29	98 nt before exon 29	ACCGTGGAGGTC-TGGAAACTCTGG	0.8 kb
418	15459-L28687	Exon 31	11 nt before exon 31	ACGCTTGTCACT-TCGGCTTCTAGG	0.5 kb
150	20227-L28154	Exon 33	2273-2274	TTCCCAGGTGAA-CGTGGCTCTCCC	0.4 kb



Length	SALSA MLPA	COL2A1	Ligation site	Partial sequence ^b (24 nt	Distance to
(nt)	probe	exon ^a	NM_001844.5	adjacent to ligation site)	next probe
337 Ж	20239-SP0944- L27954	Exon 34	2418-2419 and 2441-2442	TGCCTGGCGAGA-23 nt spanning oligo-CCCAAAGGCGAC	0.4 kb
289	15289-L27875	Exon 35	35 nt after exon 35	CTGCTGGGCATT-AGGATCCTAGCC	0.7 kb
273 Δ	20235-L27566	Exon 37	2568-2569	ACCCTCAGGGAG-AAGTTGGACCTC	0.7 kb
353	07402-L27881	Exon 39	27 nt before exon 39	CCTGCCCCTCAT-TCACCTGCTTCC	0.6 kb
130	20226-L27557	Exon 40	2796-2797	CTACTGGAGTGA-CTGGTCCTAAAG	0.5 kb
361 Ж	20242-SP0947- L27882	Exon 41	30 nt before exon 41 and intron 40-2835	TGAGGGCTTGAG-30 nt spanning oligo-GGAGCCACTGGA	0.8 kb
409	20245-L27576	Exon 42	2888 - 2889	TTCTCCTTCTAG-GGCAACCCTGGA	0.4 kb
445	20247-L27964	Exon 43	3090-3091	GTCTGGCTGGTC-AGAGAGGCATCG	0.3 kb
395 Ж	20243-SP0948- L28915	Exon 44	3192-3193 and 3227-3228	CTCCTGGAGCAT-35 nt spanning oligo-CCTGGCCTGACG	0.6 kb
213	07404-L27958	Exon 46	34 nt before exon 46	TCTTCTGGAACA-TTCTTCTCTGAG	0.7 kb
346 Ж	20241-SP0946- L27572	Exon 48	3587-3588 and 26 nt after exon 48	CTGCCCGGCCCT-29 nt spanning oligo-TCCCGAGGCCTC	0.2 kb
225	07405-L28736	Exon 49	3597-3598	CACAGGGTCCTT-CTGGAGACCAAG	0.5 kb
305 Ж	20463-SP0949- L28266	Exon 50	57 nt and 27 nt before exon 50	CCAGCCCCAGCG-30 nt spanning oligo-TTTCCTCACGAC	0.1 kb
494	20253-L27584	Exon 50	2 nt after exon 50	CGGCCCTGCTGT-AAGTGTCCTGAC	0.6 kb
238	07406-L27872	Exon 51	4007-4008	CGCACCTGCAGA-GACCTGAAACTC	0.5 kb
374	15456-L27884	Exon 52	11 nt before exon 52	GGTGCTGCCTCT-TCCCCCTGCAGG	0.4 kb
312	15455-L27878	Exon 53	99 nt before exon 53	CCCCAGTACCCT-TGAGGTCCTTGA	1.1 kb
325	07407-L28267	Exon 54	4660-4661	TTGCAAACCCAA-AGGACCCAAGTA	0.5 kb
427	15460-L28026	Exon 54	67 nt after exon 54	AGGCATGCCCAA-ATAGCAGTCCTA	
		stop codon	4617-4619 (Exon 54)		

^a See section Exon numbering on page 1 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

 Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

X This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Related SALSA MLPA probemixes

P381/P382 COL11A1 Probes for *COL11A1*, involved in type 2 Stickler syndrome and Marshall syndrome.

References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligationdependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Tarpey PS et al. (2013). Frequent mutation of the major cartilage collagen gene COL2A1 in chondrosarcoma. *Nat Genet.* 45:923-926.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P214 COL2A1

- Čopíková J et al. (2020). Expanding the phenotype spectrum associated with pathogenic variants in the COL2A1 and COL11A1 genes. *Ann Hum Genet*, 84(5), 380-392.
- Hoornaert KP et al. (2010). Stickler syndrome caused by COL2A1 mutations: genotype-phenotype correlation in a series of 100 patients. *Eur J Hum Genet*. 18:872-880.
- Kritioti E et al. (2020). First reported case of Steel syndrome in the European population: a novel homozygous mutation in COL27A1 and review of the literature. *Eur J Med Genet*, 63(7), 103939.
- Martin H et al. (2022). From First to Second: How Stickler's Diagnostic Genetics Has Evolved to Match Sequencing Technologies. *Genes*, 13(7), 1123.
- Nagendran S et al. (2012). Somatic mosaicism and the phenotypic expression of COL2A1 mutations. *Am J Med Genet A*. 158A:1204-1207.
- Richards AJ et al. (2010). Stickler syndrome and the vitreous phenotype: mutations in COL2A1 and COL11A1. *Hum Mutat*. 31:E1461–E1471.
- Wang X et al. (2016). Mutation survey and genotype-phenotype analysis of COL2A1 and COL11A1 genes in 16 Chinese patients with Stickler syndrome. *Mol Vis.* 22:697-704.

P214 prod	P214 product history		
Version	Modification		
C2	One reference probe has been replaced, one reference probe has been removed, and the sequence of one probe has been adjusted.		
C1	A total of 23 new target probes have been included and five target probes have been removed. Furthermore, all reference probes have been replaced.		
B2	One reference probe has been removed and three reference probes have been replaced. In addition, the control fragments have been adjusted (QDX2).		
B1	Twelve <i>COI2A1</i> probes, seven reference probes, and two control fragments at 100 and 105 nt have been added. Six reference probes have been replaced.		
A1	First release.		

Implemented changes in the product description

Version C2-01 - 10 August 2023 (04P)

- Product description rewritten and adapted to a new template.
- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).

Version C1-01 – 21 November 2019 (02P)

- Product description rewritten and adapted to a new template.
- Chromosomal position of the 500 nt reference probe 14894-L27890 corrected in Table 1.
- Ligation sites of the probes targeting the COL2A1 gene updated according to new version of the NM_ reference sequence.
- Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products.

More inform	More information: www.mrcholland.com; www.mrcholland.eu		
	MRC Holland bv; Willem Schoutenstraat 1		
	1057 DL, Amsterdam, The Netherlands		
E-mail	info@mrcholland.com (information & technical questions)		
	order@mrcholland.com (orders)		
Phone	+31 888 657 200		

SALSA[®] MI PA[®]